



# RibLoc® U+ 90

## Instrumentation And W&H System



### General Notes:

- The use of protective clothing (tear-resistant gloves, goggles and mask) is recommended.
- Comply with detailed information in the manufacturer's operating instructions.
- The disinfectant used must not contain chlorine.
- Use clean, dry air.
- W&H recommends the use of demineralized water.

**Review W&H IFUs for reprocessing instructions of the contra-angle handpiece, and motor with cable should always be reviewed before cleaning processing.**



### Step 1: Disassemble handpieces

- Hold contra-angle head firmly with one hand and twist sheath to release and remove.
- Remove the knee from the contra-angle sheath.
- Remove the middle gear from the knee.



### Step 2: Clean with brush and rinse with water – Manual Cleaning

- Rinse and brush off under demineralized water.
- Remove any liquid residues (absorbent cloth, blow dry with compressed air).

**⚠ Do not place the contra-angle handpieces in liquid disinfectant or in an ultrasonic bath.**

### Step 3: Reassemble Handpieces

- Insert middle gear into the knee.
- Insert the knee into the contra-angle sheath.
- Insert the sheath into the contra-angle head.
  - Note the symbols and turn until locked (see images).



### Notes:

- Types and Serial Numbers must be identical (inside contra-angle and outside on the sheath).
- Machine wash the reassembled handpieces following the manual cleaning.
- Load into W&H Tray.

### Step 4: Machine Clean - SPD

- See ACUTE Innovations IFU for details.

### Step 5: Lubricate Handpieces

- Shake and spray lubricant can before each application.
- Insert the spray nozzle into the bottom of the contra-angle handpiece.
- Cover head of handpiece with tissue, hold firmly and spray through for one second.

### Step 6: Autoclave

- See ACUTE Innovations IFU for details.